RDC-159 SOFTENING POINT METTLER

GRANULAR Material type FIECTRODES LINING



The softening point of the binder corresponds to its change from solid to liquid state. It is one of the main properties to qualify the binder, as it is usually a production target for the pitch producer and impacts most of the other binder properties. Its level should be chosen to remain below the limits of the thermal capacity of the paste plant to ensure optimum processing conditions.

The measurement is conducted with the RDC-159 apparatus, where a pressed sample (prepared with the RDC-182 apparatus) is placed in a preheated furnace. The temperature is increased with a given heat-up rate until the softened material interrupts a light barrier below the sample. The softening point Mettler is directly read from the display of the equipment and reported in °C.

| Technical information | Standard Method: | 150 5940-2 |
|-----------------------|--------------------------------------|------------------------------------|
| | Property: Softening Point Mettler | [°C] |
| | Maximum Temperature: | 400°C |
| | Sample: | Pitch pressed in a cup |
| | Process Time: | ~ 15 minutes |
| | Installation: | Workbench |
| | Dimensions (LxWxH): | 20 x 32 x 22 cm |
| | Weight: | 5 kg |
| | Electrical Property: | 230 V 1/N/PE, 50 Hz 0.9 kW, 4 A |
| | Certified Reference Material: | RDC 1159 |
| | Database Connection: | No |

Additional Recommended Equipment:

Mettler Sample Press (RDC-182)



RDC 1159



Mettler Toledo TM