

## RDC-182 Mettler Sample Press

This apparatus is used to produce the compressed sample under pre-determined conditions which is used in the softening point determination of pitch (RDC-159 Softening Point Mettler)

The pressing of the pitch into the RDC-159 cup allows the highest possible analysis precision figures as there is no bias related to melting.



\*Images and illustrations are not contractual.

<b>Specifications</b>	Measurement	Sample Preparation for RDC-159
	Sample	Fraction 2 - 1 mm
	Sample / test	1
	Process time	~ 5 minutes
<b>Configuration</b>	Set up	Workbench
	Dimensions	40 x 40 x 44 cm (LxWxH)
	Weight	35 kg
<b>Facilities</b>	Electrical connection	230V 1/N/PE, 50/60Hz
	Power	0.40 kW